



SuperEar®

**Assistive Listening Device
Model SE4000, SE5000,
SE7500 and SE9000**



Manufacturer's Suggested Cleaning Recommendations

The SE4000, SE5000 and SE7500 is packaged with Stereo Headphones with one pair of foam ear pads in place on the headphones.

Optional Headphones are available for separate purchase with the SE9000. Additional quantities of foam ear pads may be purchased from the manufacturer for \$12.95 / 20 pair. Optional/Replacement Headphones are also available for \$12.95. For infection control, sanitary Ear Pad covers are available in bulk packs of 100 covers for \$20.00.



**SE4000 SuperEar
(without Windsock)**

In a hospital or long-term care environment, the foam windscreen covering the boom microphone (SE4000 only) may be discarded, if desired. This will have no impact on sound quality, but will permit easier cleaning. After each use by a patient, the Ear Pad Covers can be removed and discarded. If desired, the headphones and transmitter may be cleaned with a light application of an alcohol-based disinfectant before replacing with new Ear Pad Covers.



**Ear Pad
Covers**

It is always best to consult with your Infection Control Officer regarding proper cleaning procedures.



SonicTechnology®

P.O. Box 539, Grass Valley, CA 95945
(530) 272-4607 □ (800) 247-5548 □ FAX (530) 272-4257
www.sonictechnology.com