



Deluxe Connected Blood Pressure Monitor UA-651BLE



Introducing the Deluxe Blood Pressure monitor UA-651BLE, one of the most advanced monitors available designed for ease of use and accuracy with a simple, compact and modern design. The UA-651BLE offers advanced features including Pressure Rating Indicator™, irregular heartbeat detection, low battery indicator, measurement time and date stamp and up to 30 readings in internal memory storage.

Unlock the Potential of the UA-651BLE with the A&D WellnessConnected™ Platform

You can now achieve greater insight into your health and wellness like never before. As part of A&D Medical's WellnessConnected™ portfolio, the UA-651BLE easily, reliably and securely connects to your iPhone or Android supported mobile device with the A&D Connect app via Bluetooth® Smart™ technology. A&D Connect is a simple, easy to use mobile application that allows you to view your wellness metrics in a whole new way, complete with trends and graphs for all of your A&D WellnessConnected™ devices, including blood pressure, weight, activity and sleep. With an A&D WellnessConnected™ account, you have unlimited storage in the A&D Cloud.

Track your wellness. Share with others.

A&D Connect allows users to engage physicians and caregivers in their health and wellness monitoring. Set your goals, see your trends and live well. The A&D Connect sharing function allows you to share your health and wellness data with family, friends and caregivers via email, social media platforms and more.



A&D Connect is compatible with Bluetooth® Smart Ready™ mobile devices from Apple and many Android device manufacturers.

APPLE:

iOS version 6 or newer




iPhone 4s or newer, iPad 3rd & 4th generation or newer, iPad Mini

ANDROID:

OS version 4.4 or newer

Requires Bluetooth Smart Ready™ device (Bluetooth 4.0).

For a list of supported devices, visit www.wellnessconnected.com/mobile





Deluxe Connected Blood Pressure Monitor

UA-651BLE

- Clinically Validated — using latest protocol ESH 2010
- SlimFit™ Cuff System w/ Gentle inflation — Flexible inflation gives less chance of re-inflation
- Advanced Measurement and Report
 - 30 reading internal memory storage / Unlimited reading storage with A&D Cloud account
 - Stores memory with time and date stamp
 - Multiple users via A&D Connect app
 - A&D Connect app allows for easy sharing of readings via mobile social applications, email and more.
- Risk Factor Indicators
 - Gives a snap shot of your BP classification in accordance with US JNC7 guidelines.
 - Irregular Heartbeat Feature – Alerts the user to the presence of an irregular heartbeat risk

Part Number	UPC Number	Cuff Size	<i>*Included. Note: Model UA-651BLE is not designed to use a small cuff.</i>
UA-651BLE	093764603327		
Medium SlimFit™ - UA-290*	093764602986	9.0" – 14.6" (23cm – 37cm)	
Large SlimFit™ - UA-291	093764602993	12.2" – 17.7" (31cm – 45cm)	

Included Accessories:

- Four AA batteries
- Bilingual Quick Start Card
- Bilingual Instruction Manual (English and Spanish)
- Carrying case

TECHNICAL SPECIFICATIONS

Measurement method	Oscillometric measurement		
Measurement range	Pressure:	20 - 280 mmHg	
	Systolic pressure:	60 - 279 mmHg	
	Diastolic pressure:	40 - 200 mmHg	
	Pulse:	40 - 180 beats /min	
Measurement accuracy	Pressure:	±3 mmHg	
	Pulse:	±5 %	
Power supply	4 x 1.5V batteries (R6P, LR6 or AA) or AC adapter (TB:233 - sold separately)		
Number of measurements	Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries) With pressure value 180 mmHg, room temperature 23 °C.		
Classification	Internally powered ME equipment (Supplied by batteries) / Class II (Supplied by adapter) Continuous operation mode		
Clinical test	According to ANSI / AAMI SP-10 1992		
Wireless communication	VZ (MURATA Manufacturing Co.Ltd.) Bluetooth Ver. 4.0LE BLP		

PACKAGING AND SHIPPING SPECIFICATIONS

Individual Unit (Box)	W x D x H (inches)	WT (lb.)	Individual Item Cube (cubic feet)
UA-651 BLE	5.12 x 4 x 7.64	1.52	.088
Unit Case	W x D x H (inches)	WT (lb.)	Case Cube (cubic feet)
10 units	21.65 x 11 x 8.46	15.87	1.18

WellnessConnected™ and Slimfit™ are trademarks of A&D Company Limited or A&D Medical. Apple name and logo are registered trademarks of Apple, Inc. Google, Android, Google Logo and Google Play are registered trademarks of Google, Inc. Facebook name and logo are trademarks of Facebook Inc. Twitter name and logo are trademarks of Twitter Inc. Microsoft Excel icon is a registered logo of Microsoft Corporation in the United States and/or other countries. Bluetooth®, Bluetooth Smart™ marks and logos are trademarks of Bluetooth SIG, Inc.

www.wellnessconnected.com



1756 Automation Parkway
San Jose, CA 95131
Tel. 1-800-726-3364
info@wellnessconnected.com